

# Meditherapy®

Accessories  
For Endoscopy



# CAN1 SERIES

## ERCP CANNULAS

DESIGNED AND PRODUCED BY



The CAN1 Series of ERCP cannulas offer excellent cannulation capabilities covering most clinical situations. All models feature a fluoroscopic tip providing X-ray visibility and are 0.035" guidewire compatible.



### FEATURES

- Single lumen design accepts a 0.035" guide wire for ease of use
- The pre-curved tip design allows the tip to be curved towards the common bile duct to facilitate cannulation
- Radiopaque tip marker aids easier identification of tip placement
- Visual distal markers aid awareness of depth of cannulation



STERILE

### SINGLE USE ERCP CANNULA

Class IIa CE0123

Part number	Guide wire compatibility	Tip size	Type	Pack size
CAN1-A3-18-200-35	0.035"	5Fr	Standard	10
CAN1-A5-18-200-35	0.035"	5.5Fr	Standard metal	10
CAN1-B3-18-200-35	0.035"	4.3Fr	Tapered	10
CAN1-D5-18-200-35	0.035"	5.5Fr	Bottle tip	10