



# CLIPMASTER3

## THREE-CLIP APPLICATION SYSTEM

DESIGNED AND PRODUCED BY  
**medw<sup>o</sup>r<sup>k</sup>**®

The revolutionary new ClipMaster3 makes it possible to endoscopically apply up to three clips in a single process, without having to be removed from the endoscope. This advance in technology ensures that time is saved to complete the procedure as well as providing a simpler technique for the application of multiple clips.

### FEATURES

- Three clips are arranged in a row inside the distal end of the sheath, allowing fast and accurate placement of multiple clips. There is no need to withdraw the device to change applicators or to re-load clips
- The compact design allows the use of standard endoscopes with an instrument channel of 2.8mm or larger
- The ClipMaster3 handle is ergonomic and comfortable to operate. Clip application is easily controlled with safe one-handed operation
- Strong stainless steel clip jaws provide a firm grasp of the mucosa
- Each Clip opens to 12 mm / 110° providing accurate placement onto the target site
- Supplied single-use to prevent cross-contamination
- Traceability : 3 detachable stickers accompany each ClipMaster3



STERILE

#### SINGLE USE MULTIPLE CLIP DEVICE

Class IIa CE0459

Part number	Compatible channel (mm)	Length (cm)	Description	Pack size
SPE1-X3-26-220	2.8 or greater	230	ClipMaster3: Disposable multiple 3-clip Clip fixing device for 2.8mm channel GI endoscopes or larger.	2