

www.meditherapy.com

Meditherapy®

Accessories
For Endoscopy



DIGNITY SHORTS AND PATIENT BIB

DIGNITY SHORTS

DISPOSABLE SHORTS FOR COLONOSCOPY

DESIGNED AND PRODUCED BY
medwork®

Patient discomfort and lack of dignity during colonoscopy can often be distressing. Dignity Shorts are an elegant response which take into account and respect patient's modesty.

FEATURES

- Adjustable access with two Velcro pads allowing the cover to be kept in the open or closed position
- Unisex design suitable for women and men of all sizes
- Single-use to reduce cross contamination



SINGLE USE DIGNITY SHORTS FOR COLONOSCOPY

Class I CE0123

Part number	Description	Leg length/width (cm)	Pack size
MIS1-Colo -Pants	Universal large size for patients undergoing colonoscopy	60/30	100
MIS1-Colo -Pants XXL	Extra large size for patients undergoing colonoscopy	60/40	100


PATIENT BIB

PATIENT PROTECTOR FOR GASTROSCOPY AND BRONCHOSCOPY

Single-use patient bib protects patient from stains and deposits during endoscopic procedures.

FEATURES

- Absorbent outer and waterproof inner layers protect the patient's clothing
- Flexible wire keeps the pouch open during the procedure
- Velcro neck strap adjusts to suit all patients



SINGLE USE PROTECTOR FOR GASTROSCOPY/BRONCHOSCOPY

Class CE0123

Part number	Length (cm)	Width (cm)	Pack size
MIS1-Gastro-Bib	55	60	100