

www.meditherapy.com

Meditherapy®

Accessories
For Endoscopy



SINGLE-USE EMR 360° SPRAY CATHETER




SPE1 SERIES

SINGLE-USE EMR 360° SPRAY CATHETER

DESIGNED AND PRODUCED BY
medw^ork®

The SPE1 EMR spray catheter is precision designed and engineered to deliver a powerful 360° spray mist pattern. Chromo-endoscopy can be performed in a uniform and clinically useful manner without any pooling of the endoscopic dye.



FEATURES

- Accurate 360° forward spray action for even and precise dyeing of the mucosa
- Available in a 220cm length for both gastroscopy and colonoscopy
- Luer lock fitting for easy connection
- Single use to reduce the potential for cross-contamination
- Traceability : 3 detachable stickers accompany each spray catheter



STERILE

SINGLE USE 360° SPRAY CATHETER

Class IIa CE123

Part number	Compatible channel size	Spray Type	Pack size
SPE1-C6-18-220	2.8mm channel or larger	360° forward	10