

www.meditherapy.com

Meditherapy®

Accessories
For Endoscopy



ENDOBITE HOT BIOPSY FORCEPS

ENDO Bite HOT

SINGLE-USE HOT BIOPSY FORCEPS

DESIGNED AND PRODUCED BY
medw^{ork}

By utilising the same pantograph hinge mechanism as reusable 'Hot' biopsy forceps to ensure the opening of the cups and optimal bite force, the EndoBite Hot biopsy forceps guarantee an optimal bite even when the endoscope is tightly angulated.

FEATURES

- 100% electrical insulation to reduce the risk of injury with current leakage through insulation pinholes
- Pantograph opening of the cups, with high precision mechanism traditionally used for reusable forceps, provides added bite force for increased biopsy quality
- Ergonomically designed handle provides a comfortable, easy to grip design, tactile control with enhanced security
- Sharp edged cups enable a 'cut' instead of a 'tear' biopsy
- 'Channel Guard' low friction sheath is a unique Fluoropolymer coating offering a very low coefficient of friction, facilitating easy passage of the biopsy forceps within the instrument channel. The lower friction sheath may also reduce wear to the biopsy channel wall
- Fenestrated cups for an increased biopsy volume
- Colour coded to differentiate between different forceps models
- Traceability : 3 detachable stickers accompany each pair of forceps



STERILE

SINGLE USE HOT BIOPSY FORCEPS

Class IIa CE0123

Part number	Compatible channel	Length (cm)	Spike	Colour code	Pack size
HOT1-C1-23-180-OL	2.0 – 2.8	180	No	Green	5
HOT1-C1-23-230-OL	2.8 – 4.2	230	No	Red	5