

Meditherapy[®]

Accessories
For Endoscopy



ENDOTRIPTOR

ERCP STONE CRUSHING LITHOTRIPTOR

The EndoTriptor contains a basket within an inner Teflon cannulating sheath and an outer crushing coil sheath. It is supplied sterile and is pre-assembled, eliminating cumbersome assembly. The EndoTriptor handle is autoclavable and incorporates a ratchet mechanism facilitating a smooth crushing action.

FEATURES

- Incorporates an inner Teflon sheath which can be extended 15cm beyond the outer coil sheath, facilitating cannulation and positioning of the crushing basket
- Injection port allows for the injection of contrast medium to confirm correct location
- Superior 16 sided basket design constructed of 4 braided wires enables the EndoTriptor to withstand a force of more than 50kg. A "rounder" shape of basket allows capture of larger stones and facilitates basket opening within the duct
- In the unlikely event of a failure due to the strain placed on the components during lithotripsy, specific "weak" points will fail, all of which are accessible and not patient threatening
- Basket and sheaths supplied sterile and pre-assembled, eliminating cumbersome assembly
- The EndoTriptor handle has been developed to be autoclavable, immersible and reusable. The top ratchet mechanism allows the handle to be turned slowly and without keeping hold, thus facilitating a smooth crushing action

2

STERILE

ENDOTRIPTOR ERCP STONE CRUSHER

Class IIa CE0459

Part number	Channel size	Type	Pack size
SEL-Sheath	4.2mm	Slide type with inner teflon catheter	1
SEL-H	-	Endotriptor handle with ratchet	1