

Meditherapy®

Accessories
For Endoscopy



IMP1 SERIES

OE SOPHAGEAL METAL SELF-EXPANDING STENTS

DESIGNED AND PRODUCED BY



The outstanding IMP nitinol metal stent range sets new standards in stent technology. The flexible, self-expanding stents have exceptional radiological visualisation characteristics, aiding in the precise positioning of the stents. In addition, repositioning threads allows for stent repositioning or complete removal.

FEATURES

- Nitinol construction ensures that the maximum diameter and radial force are reached at body temperature
- Superb expansion capabilities provide excellent opening ability even through difficult or torturous strictures
- Atraumatic edges are rounded to minimise possible mucosal damage
- IMP oesophageal stents have gold markers on both repositioning threads and there are platinum–iridium alloy markers in the middle and at both ends of the stent to ensure accurate placement
- Exceptional flexibility allows the stent to adapt to different anatomical conditions
- Full and partially covered versions coupled together with a range of diameters and lengths ensures an IMP stent is available for the full range of clinical situations
- ‘Tulip’ ends of the IMP oesophageal stents are slightly flared to aid anchoring the stent position and therefore resist migration. In addition the proximal ‘tulip’ will form a close seal with the oesophagus and therefore help ‘guide’ food through the stent
- Two radiopaque external makers are supplied with each stent



STERILE

IMP1 OE SOPHAGEAL EXPANDABLE STENTS

Class IIa CE0459

Part number	Description	Pack size
IMP1-01-20-24-080	Partly Coated (50mm). Diameter 20mm, Tulip 24mm, Length 80mm	1
IMP1-01-20-24-100	Partly Coated (70mm). Diameter 20mm, Tulip 24mm, Length 100mm	1
IMP1-01-20-24-120	Partly Coated (90mm). Diameter 20mm, Tulip 24mm, Length 120mm	1
IMP1-01-20-24-140	Partly Coated (110mm). Diameter 20mm, Tulip 24mm, Length 140mm	1
IMP1-01-24-28-080	Partly Coated (50mm). Diameter 24mm, Tulip 28mm, Length 80mm	1
IMP1-01-24-28-100	Partly Coated (70mm). Diameter 24mm, Tulip 28mm, Length 100mm	1
IMP1-01-24-28-120	Partly Coated (90mm). Diameter 24mm, Tulip 28mm, Length 120mm	1
IMP1-01-24-28-140	Partly Coated (110mm). Diameter 24mm, Tulip 28mm, Length 140mm	1
IMP1-03-24-28-060	Completely Coated. Diameter 24mm, Tulip 24mm, Length 60mm	1
IMP1-03-24-28-080	Completely Coated. Diameter 24mm, Tulip 24mm, Length 80mm	1
IMP1-03-24-28-100	Completely Coated. Diameter 24mm, Tulip 24mm, Length 100mm	1
IMP1-03-24-28-120	Completely Coated. Diameter 24mm, Tulip 24mm, Length 120mm	1