

0120912
1234546521.576232123223152
00245502687
001
2146
0142
1024

Meditherapy®

Accessories
For Endoscopy



12301012480701224655-57666
01547893468754642213520216
14441
016623284985322.20353548
0123412
1234546521.5762321
1234546521.5762321
00245502687
001
2146
3200
01224555457886015478934

1234546521.5762321

INJ1 SERIES INJECTION NEEDLES

INJ1 SERIES

INJECTION NEEDLES

DESIGNED AND PRODUCED BY

medwork®

The INJ1 injection needle range incorporates a 'Lanzetten' cut, which is an extra sharp triple bevel ensuring optimal puncture every time. The handle's double 'click-stop' mechanism, allows the needle to be locked in the extended position or be fixed in the retracted position minimising the possibility of needle-stick injury.



FEATURES

- Lanzetten cut triple-bevelled needle for a highly effective puncture
- Ergonomic handle with a reversible 'click / unlock' system which maintains the position of the needle in either the extended or retracted positions providing extra safety and scope protection
- Transparent distal sheath helps visualise the needle position
- Easy to operate with one hand
- Distal sheath reinforced with a radio-opaque metallic tip providing excellent puncture ability



STERILE

INJ1 SERIES – SINGLE USE INJECTION NEEDLE WITH METALLIC TIP

Class IIa CE0123

Part number	Channel size	Working length (cm)	Gauge / Ø needle (mm)	Needle length (mm)	Pack size
INJ1-A1-05-5-18-180	2.0mm or larger	180	25g/0.5	5	10
INJ1-A1-07-5-18-180	2.0mm or larger	180	25g/0.5	5	10
INJ1-A1-07-5-23-180	2.8mm or larger	180	23g/0.7	5	10
INJ1-A1-07-5-23-230	2.8mm or larger	230	23g/0.7	5	10
INJ1-A1-09-5-23-180	2.8mm or larger	180	21g/0.9	5	10
INJ1-A1-09-5-23-230	2.8mm or larger	230	21g/0.9	5	10



STERILE

INJ1 SERIES – SINGLE USE FIBRIN/THROMBIN DUAL NEEDLE WITH METAL TIP

Class IIa CE0123

Part number	Channel size	Working length (cm)	Gauge / Ø needle (mm)	Needle length (mm)	Pack size
INJ1-D1-08-5-26-210	2.8mm or larger	210	21g/1.0	5	1



STERILE

UNO-JECT- S – SINGLE USE EMR INJECTION NEEDLE WITH TAPERED TIP

Class IIa CE0123

Part number	Channel size	Working length (cm)	Gauge / Ø needle (mm)	Needle length (mm)	Pack size
AS2050523	2.8mm or larger	230	25g/0.5	5	5