

# Meditherapy<sup>®</sup>

Accessories  
For Endoscopy



# PRO1 SERIES

## PLASTIC BILIARY STENTS

DESIGNED AND PRODUCED BY



The PRO1 Series incorporates a wide range of biliary plastic stents offering various shapes and size options. Their optimal combination of different sizes and designs allows for the most appropriate selection of stent for every eventuality.



### FEATURES

- Manufactured from ultra-smooth polyethylene material for enhanced cannulation features
- Flexible design helps reach hilar strictures and successfully deploy at difficult angles
- Fluoroscopic visual marker bands to aid correct placement
- Full range including pre-assembled individual pushers and guiding catheters for every possibility



STERILE

### PRO1 SERIES BILIARY PLASTIC STENTS - (Most common specifications) Class IIa CE0123

Part number	Diameter	Length	Type	Pack size
PRO1-A3-10-04-06	10Fr	4cm	Straight	1
PRO1-A3-10-06-08	10Fr	6cm	Straight	1
PRO1-A3-10-09-11	10Fr	9cm	Straight	1
PRO1-A3-10-12-14	10Fr	12cm	Straight	1
PRO1-D4-08-06-09	8.5Fr	6cm	Pigtail	1
PRO1-D4-08-09-12	8.5Fr	9cm	Pigtail	1
PRO1-D4-08-12-15	8.5Fr	12cm	Pigtail	1
PRO2-F5-05-220-25	8.5Fr guide catheter	220cm	Guide catheter for 8.5Fr stents	1
PRO2-G7-08-170	8.5Fr push catheter	170cm	Push catheter for 8.5Fr stents	1
PRO2-F5-06-220-35	10Fr guide catheter	220cm	Guide catheter for 10Fr stents	1
PRO2-G7-10-170	10Fr push catheter	170cm	Push catheter for 10Fr stents	1
PRO1-H7-05-205-25	For 8.5Fr stents	220cm	8.5Fr pre-assembled introduction kit	1
PRO1-H7-06-205-25	For 10Fr stents	220cm	10Fr pre-assembled introduction kit	1

Full range of French sizes and stent shapes available upon request

Birchwood House • 55 Vanguard Way • Shoeburyness • Essex S53 9QY • UK

Tel: +44 (0)1702 291878 • Fax: +44 (0)1702 290013 • Email: info@medinn.co.uk • Web: www.meditherapy.com