

# Meditherapy®

Accessories  
For Endoscopy




UNO-DRAIN ENBD TUBES

# UNO-DRAIN

## NASAL BILIARY DRAINAGE TUBE

The Uno-Drain ENBD tubes provide drainage of biliary ducts and/or the injection of solutions. Uno-Drain is also available in kits, grouping all items needed to provide easy and rapid placement.



### FEATURES

- Axial and lateral holes in distal end provide easy diffusion or suction of liquid
- Stainless steel distal ring for excellent radiopacity
- Made from specialist materials which give tube excellent memory shape and make it easy to pass over a guidewire

### KIT COMPOSITION:

- 1 naso-biliary tube
- 1 nasal tube Ø 5, 7 or 10 Fr / length 50 cm
- 1 extension containing a 3-way stopcock, length 135 cm
- 1 removable luer-lock connector Ø 5, 7 or 10 Fr
- 1 guidewire Ø 0,035" or 0,025" / length 400 cm



STERILE

### SINGLE USE NASAL BILIARY DRAINAGE TUBE

Class IIa CE0459

Part number	Curve type	Ø Tube (Fr)	Length (cm)	Reference of corresponding kit	Pack size
AS0150529	Curve 1 : short straight tapered loop	5	290	AS0460500	1
AS0150729	Curve 1 : short straight tapered loop	7	290	AS0460700	1
AS0091029	Curve 1 : short straight tapered loop	10	290	AS0461000	1
AS0160529	Curve 2 : double loop pigtail	5	290	AS0470500	1
AS0160729	Curve 2 : double loop pigtail	7	290	AS0470700	1
AS0101029	Curve 2 : double loop pigtail	10	290	AS0471000	1