

# Meditherapy®

Accessories  
For Endoscopy




UNO-TRACT BILARY BALLOONS

# UNO-TRACT

## BALLOON EXTRACTION CATHETER

By combining innovative design features and high quality materials, the Uno-Tract is one of the best balloon extraction catheters available.



### FEATURES

- Excellent balloon strength to guarantee efficiency and security of extraction
- Distal metal markers for fluoroscopic identification of balloon position
- Coloured graduation markings on distal end give excellent visualisation during video-endoscopy procedures
- 18mm version available for dilated ducts
- Tapered model 7 Fr to 5 Fr for ease of cannulation
- Triple lumen model: 1 for guidewire / 1 for inflating-deflating of the balloon / 1 for injection



STERILE

### SINGLE USE EXTRACTION BALLOON CATHETER 2 LUMEN

Class IIa CE0459

Part number	Ø catheter (Fr)	Length (cm)	Ø inflated balloon (mm)	Recommended Guidewire	Pack size
AS0030520	5	200	12	0,025"	5
AS1910720	7	200	12	0,035"	5
AS0030720	7	200	16	0,035"	5
AS0280720	7	200	18	0,035"	5
AS1100720	7 => 5	200	16	0,035"	5

### 3 LUMEN

AS1420720	7	200	16	0,035"	5
-----------	---	-----	----	--------	---