

Meditherapy[®]

Accessories
For Endoscopy



UNO-TRAP POLYP RETRIEVAL DEVICE


UNO-TRAP

FILTERED CONTAINER FOR POLYP RETRIEVAL

The Uno-Trap helps protect staff by minimising contact with resected polyp specimens while conserving sample integrity. The Uno-Trap is an essential device for polypectomy procedures.

The Uno-Trap comprises:

- Container divided into 4 identical compartments labelled 1 to 4
- Lid containing 2 tubes, a central tube with suction attachment and another to attach to endoscope



FEATURES

- Single-use device which captures resected polyps in closed system to avoid risks of cross infection and contamination
- Classification of 4 compartments to preserve chronology of resections
- Fluorescent marking enables easy identification of compartment numbers, even in darkened endoscopy rooms
- Lid is sealed with one tube to ensure retrieved polyps are transported to laboratory without risk of samples being damaged



NON STERILE

SINGLE USE CONTAINER FOR POLYP RETRIEVAL

Part number

Description

Class I CE

Pack size

10017130B

Non sterile units individually packaged

20