

Meditherapy®

Accessories
For Endoscopy



UNO-VALVE

DISPOSABLE STERILE ENDOSCOPE BIOPSY VALVE

The Uno-Valve range of disposable endoscopic biopsy valves eliminates the concern of cross-contamination between patients and minimises nurse exposure to potentially infectious body fluids due to reused or worn out biopsy valves. The valves are individually sterile packed for enhanced safety.

FEATURES

- Single use to minimise the risks associated with cross-contamination
- Individually sterile packaged for enhanced safety
- Valves have fluorescent properties which enable them to be seen in darkened procedurerooms, therefore ensuring accurate placement of accessories
- Suitable for Olympus, Pentax and Fujinon endoscopes
- Close fitting cap prevents leaking of fluids
- Conforms to BSG Guidelines advising the use of disposable valves



UNO-VALVE DISPOSABLE BIOPSY VALVE			Class I CE0459
Name	Part number	Type	Pack size
Uno-Valve	AS2430000	Uno-Valve: Disposable biopsy valve for Olympus and Fujinon GI endoscopes.	100
Uno-Valve	AS2440000	Uno-Valve: Disposable biopsy valve for Pentax GI endoscopes.	100
Uno-Valve Plus	AS2450000	Uno-Valve Plus: Disposable biopsy valve with irrigation adaptor for Olympus and Fujinon GI endoscopes.	100
Uno-Valve Plus	AS2460000	Uno-Valve Plus: Disposable biopsy valve with irrigation adaptor for Pentax GI endoscopes.	100