

# Meditherapy®

Accessories  
For Endoscopy

010 1 01 1



010 1 01 1

10 01 01011 0101 01



# UNO-VIEW

## OESOPHAGEAL VARICES LIGATING DEVICE

The Meditherapy Uno-View multiple band ligator is a disposable ligation device for oesophageal varices. The ligating unit is designed to deploy six latex bands independently and successively from a pre-loaded cylinder.

### FEATURES

- Minimal strings and knots provides improved visualisation
- Non-elastic trip wire ensures the bands are released as the handle is turned on a one-to-one basis. This also does not put traction onto the scope and therefore does not cause 'tilting' of the endoscope tip
- Coloured fifth band aids identification of number of bands placed
- Only one barrel thread which provides a less fussy, clearer view
- Simple handle design with lock/unlock feature improves user operatorability
- Suitable for 9.5 to 12mm endoscopes ensuring proper fit with a wide range of scope sizes
- The slim, short barrel makes the Uno-View more patient tolerant and easier to intubate
- A flushing needle is provided to flush the biopsy channel during band ligation



UNO-VIEW OESOPHAGEAL VARICES LIGATOR			Class I CE
Part number	Description	Size of scope	Pack size
SBLS-6-0	Six shot multiple band ligator for oesophageal varices. For Olympus and Pentax scopes	9.5-12mm	1
SBLS-6-F	Six shot multiple band ligator for oesophageal varices. For Fujinon scopes	9.5-12mm	1
SBL	Spare barrel of 6 bands only	9.5-12mm	1