

Meditherapy®

Accessories
For Endoscopy



CLEANSCOPE® STORAGE & TRANSPORTATION

Available exclusively through

ASP ADVANCED STERILIZATION PRODUCTS®
a *Johnson & Johnson* company

CLEANASCOPE®

SECURE TRANSPORTATION AND SHORT TERM STORAGE OF ENDOSCOPES

Cleanascope® is a unique transport and short term storage solution for endoscopes and other delicate equipment, which minimises direct handling and airborne contamination.

FEATURES

- Safe and easy handling of endoscopes and other delicate equipment
- Hygienic short-term storage of clean scopes before endoscopy list starts
- Helps prevent airborne contamination of disinfected scopes before procedure
- Hygienic storage of scopes whilst awaiting disinfection after use
- Preventing contamination from dirty endoscopes after procedure
- Clear identification of 'clean' or 'contaminated' equipment
- Helps avoid accidental damage to your scopes thereby reducing repair bills.

CLEANASCOPE® SIX TRAY CART Class I CE

Reference	Description	Height	Width	Depth
103630	Six tray cart + six transport trays	925 mm	520 mm	570 mm

TEN TRAY CART

103631	Ten tray cart + ten transport trays	1435 mm	520 mm	570 mm
--------	-------------------------------------	---------	--------	--------

TWO TRAY TRANSPORT CART

103632	Two tray cart + two transport trays	925 mm	520 mm	570 mm
--------	-------------------------------------	--------	--------	--------

FIVE TRAY CART

103633	Five tray cart + five transport trays	800 mm	520 mm	570 mm
--------	---------------------------------------	--------	--------	--------

TRANSPORT TRAY

Reference	Internal/external Length	Internal/external width	Internal/external depth	Unit
103614	475 / 520 mm	370 / 420 mm	95 / 100 mm	1

WALL MOUNTED DISPENSER FOR DISPOSABLE LINERS

Reference	Description	Unit
103604	Wall mounted dispenser for disposable liners	1

CLEANASCOPE® : SINGLE USE LINERS



Reference	Description	Unit
103580	1 set individually packaged	200



Coronation Road, Ascot, Berkshire, SL5 9EY. **Customer Services:** 0800 864 060
Email: contact@medgb.jnj.com